

SUPER GROW PROFESSIONAL HOME LASERS

INSTRUCTIONS

Covering all Super Grow Laser models
50 lasers to 200 lasers and all others in between
Dome and Panel systems

GENERAL USE

Place chair under the stand and adjust the head so the lasers are about one or two inches from the top of the scalp. Keep the laser a half inch to two inches from scalp.

For the first week, use every day for twenty minutes. Then use for twenty minutes every other day for two months. After that you can use three times per week. See the chart.

It is normal for shedding to occur as your follicles go through a period of new growth. This will not last long.

Do not use every day after the first week. The follicles will stop responding to the laser.

It is best to use on a clean scalp. If you use topical products, use them after using the laser.

It can take up to four to six months to see very positive results. Don't get discouraged and give up.

GAPS IN USE

If you have any gaps in use, because of leaving town or some other reason, this is okay. Just resume using the laser on a regular schedule when you return. You do not need to start over from the beginning.

WARRANTY

Your Super Grow Laser is covered for one year from the date of purchase. If the item malfunctions, send us an email for shipping instructions and we will repair or replace it. After the warranty period we will fix it for you at nominal charge for repair.

Tampering with or modifying the unit in any way voids the warranty.



SUPER GROW LASERS USAGE CHART

	Sun	Mon	Tues	Wed	Thurs	Fri	Sat	
Week 1	Red	Red	Red	Red	Red	Red	Red	Every Day 1st Week
Week 2	White	Red	White	Red	White	Red	White	
Week 3	Red	White	Red	White	Red	White	Red	Every Other Day After First week
Week 4	White	Red	White	Red	White	Red	White	
Week 5	Red	White	Red	White	Red	White	Red	
Week 6	White	Red	White	Red	White	Red	White	
Week 7	Red	White	Red	White	Red	White	Red	
Week 8	White	Red	White	Red	White	Red	White	
Continue on this schedule which repeats every two weeks	Red	White	Red	White	Red	White	Red	
	White	Red	White	Red	White	Red	White	